

# CHISON

Value Beyond Imaging

## CBit 9 *PRO*

Everything you need





## 10.1 inch HD touch panel

Superbly responsive

Customizable interface

Ergonomic tilting ensures all-dimensional,  
multi-angle visualization



# Ergonomic design

21.5 inch Wide HD LED

Gel warmer


10.1 inch Touch screen

Height adjustable

Battery

Four probe connectors





## Dual-Screen Display


Dual-screen display for better view: touch screen, main screen

Easy to slide

Finger- gesture to ZOOM in and ZOOM out

Convenient measurement ,rotation

## Highly Improved productivity




### Intelligent Doppler

- Automatically adjust the ROI direction, PFR in color mode, doppler gate in PW mode.
- More efficient.
- Much easier to measure.
- Allow you to focus on patients better.



### Intelligent Focus

- Automatically optimize the focus position.
- Focus on the interesting area to improve the image.
- Efficient and intelligent.




### Auto Breast Lesion Detection

- Automatically detect breast lesion
- Provide the size
- Efficient for the diagnose




# Expands the speciality

## SonoContrast




- An innovative technology.
- Higher sensitivity.
- Better contrast.
- Adapt to all probes: linear, convex, phased array, transvaginal.

## SonoNeedle




- Enhance the display of needle during puncture.
- Better performance than traditional needle enhance technology.
- No angle dependence.

## Color Panoramic imaging



- Enlarge the scanning area in real time, provide the more color information.
- Useful for the large vessel especially the upper extremity and low extremity.


## Auto Thyroid



- Automatically locates and traces the thyroid Cyst.
- Automatic measures the size, area and perimeter.


## Raw Data

- Adjust the parameters after freeze (gain, TGC, dynamic) .
- Speedy the scanning time, save the processing time.
- Efficiency and fast .




## Intelligent-OB

- Intelligent software for OB, High efficiency and precise measurement tools.
- Automatically measure: BPD, HC, AC, FL, NT.
- One step to obtain the result.



## SonoCrystal

- To do prenatal diagnosis especially assess the spine, ribs and skin.
- Effective enough to differentiate soft tissue and bony structures.
- Give particular information on the contour of the ribs and enabled inference on mineralization.




## SonoColor


- Automatic measure the velocity of blood flow in C Mode.
- Provide quantitative flow analysis without spectrum.




# Better diagnostic confidence



Carotid, B Mode



Jugular Vein, B Mode




Kidney, C Mode




Left ventricular short axis, B Mode




Left ventricular long axis, B Mode




Thyroid, B Mode



MSK, B Mode



Vertebral Artery, B Mode



Fingertip blood flow, C Mode



Convex  
D3C60L



Linear  
D12L40L



Linear  
D7L40L



Volume  
V4C40L



Transvaginal  
D7C10L



Phased array  
D3P64LS



Phased array  
D3P64L



Phased array  
D5P64L



Trans-Rectal  
D7L40L-REC



Micro-Convex  
D6C15L

## CHISON Medical Technologies Co., Ltd.

**Sales & Service Contact Address:** No.3, Changjiang South Road, Xinwu District, Wuxi, Jiangsu, China 214028

**TEL :** 0086-510-85310937 **FAX :** 0086-510-85310726 **EMAIL :** export@chison.com.cn

We reserve the right to make changes to this catalogue without prior notice. All rights reserved.  
Please contact our local dealer for the latest information.